

Thrombocyte Shaker

Keeps your thrombocytes in motion



Thrombocyte Shaker 60A



Thrombocyte Shaker 96A



Thrombocyte Shaker 96



Thrombocyte Shaker 60



Thrombocyte Shaker 48



Thrombocyte Shaker 24

Storage of thrombocytes – standstill or motion

- Upper section moves over 30mm at ~60 strokes at 50Hz / ~72 strokes at 60Hz per minute
- Every shelf has a capacity of 6 bags (for models 24, 48, 60, 60A)
- Every shelf has a capacity of 12 bags (for model 96, 96A)
- Built in drawer stopper prevents accidental falling out
- Designed in stainless steel
- Each drawer can be operated independently from other drawers
- Designed for silent and reliable long time use
- No-motion alarm
 - Acoustic and visual warn signal
 - Connector for external no motion alarm

Technical data	Thrombocyte Shaker 24	Thrombocyte Shaker 48	Thrombocyte Shaker 60	Thrombocyte Shaker 96	Thrombocyte Shaker 60A	Thrombocyte Shaker 96A
Capacity (bags)	24	48	60	96	60	120
Capacity (Apheresis bags)					20	50
Dimensions DxWxH (cm)	38,1x44,8x25,9	38,1x44,8x38,8	38,1x44,8x45,3	33,8x87,7x39,3	37,7x47,2x45,3	37,6x99,6x45,7
Weight (kg)	20	25,5	28,25	37,2	28	53
Voltage	110-120V\220V-240V and 50HZ/60Hz					
Compliance data						
Classification	Medical device class I					
Directive and Standards applied	MDD 93/42/EEC EN 60601-1; EN 60601-1-2 ; RoHS2					
Classification (Protection against electric shock)	Class I, Type B					